

- Type I vs Type II -



<Type I> RWS/L-42, 48, 54

★ Easy fit to the dental arch

Anterior plastic part can easily fit to the dental arch with a hair dryer

★ Wide Variety


42mm, 48mm, 54mm

Choose according to the length between canines.

- Comparison Type I with Type II -
(Extracted dental model is used)



Type II RWL-P50



Type I RWL-42



Type I RWL-48



Type II RWL-P36

<Type II> RWS/L-P36, P50

★ Strengthened

Improved strength than Type I

(It is a bit difficult to fit to the dental arch in comparison with Type I.)

★ Convenient and priced low

Longer plastic part (50mm),

But same positioning of loop as RWL-42(42mm)

A more beautiful smile line for you...

→It is usable for patients

who use 42mm or 48mm size retainer.

→It is suitable for extraction cases as well.

New!

P36 is suitable for lower arch.

(Every type has looped and non-looped version)